

**One-step breast reconstruction** with  
**Pectoralis Major intact**

**BRA****ON**<sup>®</sup>  
Muscle sparing one-step reconstruction

# The simpler the better

BRAXON® is fruit of the combined experience in the tissue engineering field and clinical practice.

In 2012, the heartfelt necessity of more conservative surgical procedures boosted innovation in reconstructive surgery. Giving a nod to biomaterials with regenerative capacities, DECO med® disclosed an innovatory surgical technique to the Scientific Community and ADM-assisted one-step prepectoral breast reconstruction was born.

When indicated, prepectoral breast reconstruction provides important benefits on multiple fronts. Firstly, to the patient, who has her physical integrity recovered in a single operation while obtaining a natural cosmetic outcome. Then, to the healthcare system, by cutting the costs associated with invasive and complications-prone surgical procedures.

**BRAXON® is the all-round approach to breast reconstruction with the most advanced biomaterial and the least impact on patient's anatomy.**



# WHY

Because sparing the Pectoralis Major changes everything. Lifting it for submuscular breast reconstructions results in various functional deficits involving muscle strength and arm range of motion. Over time the interference in the synergy of the various muscle groups may result in weakness of humeral-scapular articulation with repercussions in movement, causing pain, and unsatisfying cosmetic results.

# WHAT

BRAXON<sup>®</sup> is the pre-shaped porcine dermal matrix which allows the tailoring of an ADM suit around the mammary prosthesis and its fixing on the pectoralis muscle, which is kept intact. Covered by several patents, the BRAXON<sup>®</sup> silhouette perfectly adapts to the shapes of all silicone prostheses. Publications from the last two decades have demonstrated a lower incidence of capsular contracture when the breast prosthesis is covered by a dermal matrix.

# HOW

In case of one-step/immediate breast reconstruction, the creation of a sub-muscular pocket takes longer and requires complex postoperative rehabilitation. All this is avoided with BRAXON<sup>®</sup>-assisted prepectoral breast reconstruction. Creating a complete biological envelope for the synthetic implant, BRAXON<sup>®</sup> constitutes a new interface between implant and patient's soft tissues reducing foreign-body reaction. This allows the complete subcutaneous placement of the prosthetic implant where the breast belongs: above pectoral muscles.

# WHEN

When breast reconstruction is planned following nipple- or skin-sparing mastectomy. When a good, vascularised subcutaneous layer is present, in fact, previous radiotherapy or co-existing medical conditions such as diabetes or connective tissues disorders are contraindications. It is advised to maintain a viable subcutaneous tissue during mastectomy: the use of monopolar diathermy should be reduced to a minimum to prevent skin flap necrosis.

# Natural, Artificial, Synthetic

We have been used to consider biocompatibility as a dichotomy splitting the world in what is biocompatible and what is not. Recent evidence shows how the effect of an implant varies depending on the surgical site, the aim of the surgery, the shape of the implant, and the nature of the material it is made of.

Natural matrices are made of native proteins (collagen) derived from acellularised animal tissues. The one BRAXON® is made of derives from porcine dermis: the organism recognises it as its own and transforms it into new self-tissue through the natural regenerative process (remodelling).


Artificial membranes derive from natural polymers, but they are chemically modified for reinforcement (e.g. cross-linked membranes). They can be either tolerated or digested by the organism, but they do not stimulate any regenerative process.

Synthetic meshes, made of polymers obtained by chemical synthesis (e.g. plastics like polyethylene, polylactide, polypropylene), are chemically tolerated but they do not stimulate a regenerative process, with the organism trying to confine them in a fibrous capsule.


**Nowadays the old concept of biocompatibility is not enough to define the effectiveness of a biomaterial.**

**Its performance must surpass passive tolerance.**


## Histologies



I. Haematoxylin-Eosin staining of a sterile BRAXON® sample. The staining reveals the complete absence of cellular material. 100x magnification.



II. Azan-Mallory staining of a sterile BRAXON® sample. Collagen fibres of the matrix are highlighted, absence of cellular material. 50x magnification.



III. Haematoxylin-Eosin staining of a BRAXON® sample 4 weeks after implantation. Nuclei, stained darker, show a high degree of cell colonisation. 50x magnification.

# Technical characteristics

BRAXON<sup>®</sup> is a 0.6 mm thick Acellular Dermal Matrix (\*) derived from selected porcine dermis. It is specifically designed for fast integration without amplifying the inflammatory process.

The exclusive production process has been developed with the aim of generating a completely natural product. BRAXON<sup>®</sup> is not cross-linked and is free from any chemical substance which can amplify the inflammatory response and slow the path of tissue regeneration. The native protein structure provides immediate bio-availability for incorporation into the host tissue. A fast integration lessens the inflammatory response and its effects, such as redness and seroma which commonly derive from the implantation of other biomaterials in the breast.

To preserve the genuine qualities of BRAXON<sup>®</sup> and to avoid preservatives that could amplify postoperative inflammation, in the final phases of the production a sublimation step is applied. This freeze-drying process allows the complete removal of any liquid through exact use of pressure and temperature. The result is a


reliable product, optimal for long-lasting room temperature preservation. A simple rehydration before use is all it takes, without the need for repeated washings to attempt the removal of added chemicals which other products contain.

BRAXON<sup>®</sup> patented shape allows a perfect fit around the silicone prosthesis creating a new interface that can be sutured onto the muscle surface. It fits neatly under the skin in a perfect position, which can be reliably checked prior to closing the surgical access.

After surgery, early and prolonged use of a conforming bra for 3 or 4 weeks as well as a compressive dressing and reduced mobility of the arm will significantly decrease the risk of seroma formation.


**A biomaterial must be active, not merely inert, with the aim of driving the biological process of guided tissue regeneration.**

(\*) Acellular Dermal Matrix: a complex net of extracellular macromolecules that, in addition to performing a connective function between cells and tissues, provides an organized structure in which cells can migrate and interact with each other.



**CASE OF BREAST  
IMPLANT SIZE REVISION:  
A WELL-VASCULARISED  
NEO-FASCIA WAS OBSERVED  
AFTER 7 MONTHS  
FROM THE BRAXON<sup>®</sup>  
PROCEDURE**

# Surgical Steps



**START 00'**  
BRAXON® procedure begins during mastectomy





**05'**  
BRAXON® is hydrated for 5 minutes to make it soft and pliable



**15'**  
Use of sizer to choose the right size and shape of the breast implant

## Tailoring BRAXON®

BRAXON® is pre-shaped in such a way as to contain a breast implant of any size and shape. Its use is intuitive and only requires scissors and sutures to “dress” the prosthesis and be sutured upon the pectoralis muscle.



35'

45'

FINISH

60'




**35'**  
BRAXON<sup>®</sup> tailoring around the mammary prosthesis



**45'**  
The BRAXON<sup>®</sup>-implant construct is inserted and sutured to the chest



**FINISH 60'**  
Skin is closed and an appropriate dressing is applied



# The iBAG Study

THE LARGEST MULTICENTRE DATA COLLECTION ON PREPECTORAL BREAST RECONSTRUCTION: THE iBAG STUDY

Jaume Masiá | iBAG Working Group




BRAXON® patented shape allowed for the standardised wrapping of the breast implant granting reproducible results.

This favoured the collection and analysis of 1450 cases of BRAXON® implantation from 30 international centres.

From the iBAG (international Braxon Audit Group) and with the endorsement from independent scientific organizations (INPECS and Cochrane collaboration), the widest and most solid scientific evidence on prepectoral breast reconstruction was elaborated.


The study reports low complication rates together with natural and symmetrical aesthetic results. The statistical analysis confirms the guidelines proposed for patient selection, highlighting higher risks for smokers, immunosuppressed patients, those with diabetes, and in case of previous surgery or radiotherapy. On the other hand, a new possible indication is represented by those cases with planned postoperative radiotherapy (PMRT). In an average follow-up of 2 years, BRAXON® has been proven reducing the rate of capsular contracture and improving results in terms of tissue softness.

COMPLICATIONS	NON-IRRADIATED BREASTS N =1243	BREASTS WITH PMRT N =159
Seroma	7.0%	10.7%
Dehiscence	4.5%	4.4%
Haematoma	2.2%	2.5%
Necrosis	3.1%	3.8%
Infection	4.7%	5.0%
Extrusion	1.0%	1.3%
RBS	3.5%	2.5%
Fever	1.9%	0.6%
Implant rotation	0.2%	0.0%
Capsular contracture	1.7%	6.3%
Rippling	3.1%	0.0%
Other	3.1%	4.4%
Implant loss	5.6%	10.7%



Bilateral BRAXON® implantation with anatomical implants: 15 months postoperative photographs

# THE WINNING IDEA



## References

### BRX06S

BRAXON<sup>®</sup> acellular dermal matrix.  
Pre-shaped ADM for breast implant  
total coverage.



# BRAXON® Scientific Evidence

More than 60 scientific studies demonstrate the safety and efficacy of Braxon technique exploring the different aspects of the original prepectoral surgery. Over 200 authors contributed to its scientific corpus, providing patients' selection criteria, surgical guidelines, and postoperative management tips.

1. Berna G, et al. Evaluation of a novel breast reconstruction technique using the Braxon® acellular dermal matrix: a new muscle-sparing breast reconstruction. *ANZ Journal of Surgery*. 2014. DOI: 10.1111/ans.12849
2. Maruccia M, et al. One-stage breast reconstruction techniques in elderly patients to preserve quality of life. *European Review for Medical and Pharmacological Sciences*. 2016. PMID: 28051266
3. Iqbal FM, et al. Host Integration of an Acellular Dermal Matrix: Braxon Mesh in Breast Reconstruction. *Clinical Breast Cancer*. 2016. DOI: 10.1016/j.clbc.2016.06.009
4. Vidya R, Cawthorn SJ. Muscle-Sparing ADM-Assisted Breast Reconstruction Technique Using Complete Breast Implant Coverage: A Dual-Institute UK-Based Experience. *Breast Care*. 2017. DOI: 10.1159/000464401
5. Vidya R, Iqbal FM. Breast anatomy: Time to classify the subpectoral and prepectoral spaces: Letter to the Editor. *Clinical Anatomy*. 2017. DOI: 10.1002/ca.22878
6. Berna G, Cawthorn SJ. Long term follow-up on prepectoral ADM-assisted breast reconstruction: evidences after 4 years. *European Journal of Plastic Surgery*. 2017. DOI: 10.1007/s00238-017-1285-z
7. Vidya R. Prepectoral Breast Reconstruction or Muscle-Sparing Technique with the Braxon Porcine Acellular Dermal Matrix: Plastic and Reconstructive Surgery - Global Open. 2017. DOI: 10.1097/GOX.0000000000001364
8. Linguadoca C, et al. L'HTA applicato alla ricostruzione mammaria in interventi chirurgici one-step: dall'utilizzo degli espansori all'impiego di matrici biologiche, l'esperienza dell'Azienda Ospedaliero-Universitaria di Parma. *HTA Focus - Pills of Clinical Governance*. 2017. DOI: 10.23753/htafocus2017.01.003
9. Vidya R, Iqbal FM. A Guide to Prepectoral Breast Reconstruction: A New Dimension to Implant-based Breast Reconstruction. *Clinical Breast Cancer*. 2017. DOI: 10.1016/j.clbc.2016.11.009
10. Tasoulis M, et al. Subcutaneous implant breast reconstruction: Time to reconsider? *European Journal of Surgical Oncology*. 2017. DOI: 10.1016/j.ejso.2017.04.008
11. Onesti MG, et al. Clinical, histological, and ultrasound follow-up of breast reconstruction with one-stage muscle-sparing "wrap" technique: A single-center experience. *Journal of Plastic, Reconstructive & Aesthetic Surgery*. 2017. DOI: 10.1016/j.bjps.2017.06.023
12. Vidya R, et al. Evaluation of the effectiveness of the prepectoral breast reconstruction with Braxon dermal matrix: First multicenter European report on 100 cases. *The Breast Journal*. 2017. DOI: 10.1111/tbj.12810
13. Jafferbhoy S, et al. Early multicentre experience of pre-pectoral implant based immediate breast reconstruction using Braxon®. *Gland Surgery*. 2017. DOI: 10.21037/g.s.2017.07.07
14. Maruccia M, et al. One-stage muscle-sparing breast reconstruction in elderly patients: A new tool for retaining excellent quality of life. *The Breast Journal*. 2018. DOI: 10.1111/tbj.12860
15. Vidya R, et al. Management based on grading of animation deformity following implant-based subpectoral breast reconstruction. *Archives of Plastic Surgery*. 2018. DOI: 10.5999/aps.2017.01242
16. Salgarello M, et al. Current evidences on immediate breast reconstruction after mastectomy. *Translational Cancer Research*. 2018. DOI: 10.21037/tcr.2018.03.40
17. Cattelani L, et al. One-Step Prepectoral Breast Reconstruction With Dermal Matrix–Covered Implant Compared to Submuscular Implantation: Functional and Cost Evaluation. *Clinical Breast Cancer*. 2018. DOI: 10.1016/j.clbc.2017.11.015
18. Chandarana MN, et al. Acellular dermal matrix in implant-based immediate breast reconstructions: a comparison of prepectoral and subpectoral approach. *Gland Surgery*. 2018. DOI: 10.21037/g.s.2018.03.05
19. Cattelani L, Polotto S. The Economics of Prepectoral Breast Reconstruction. *Plastic and Reconstructive Surgery*. 2018. DOI: 10.1097/PRS.0000000000004660
20. Gardani M, et al. Prepectoral breast reconstruction using the Braxon® porcine acellular dermal matrix: a retrospective study. *European Journal of Plastic Surgery*. 2018. DOI: 10.1007/s00238-018-1455-7
21. Salgarello M, et al. A case of important weight loss after a prepectoral breast reconstruction. *European Journal of Plastic Surgery*. 2018. DOI: 10.1007/s00238-018-1413-4
22. Kankam HKN, et al. Combination of acellular dermal matrix with a de-epithelialised dermal flap during skin-reducing mastectomy and immediate breast reconstruction. *annals*. 2018. DOI: 10.1308/rcsann.2018.0127
23. Vidya R, Green M. Minimal Pain with Prepectoral Implant-Based Breast Reconstruction. *Plastic and Reconstructive Surgery*. 2019;143(1):236e. DOI: 10.1097/PRS.0000000000005135
24. Chandarana M, et al. Outcomes of prepectoral implant-based breast reconstruction with Braxon® acellular dermal matrix—a single-centre experience. *European Journal of Plastic Surgery*. 2019. DOI: 10.1007/s00238-019-01512-2
25. Vidya R, et al. Prepectoral implant-based breast reconstruction: a joint consensus guide from UK, European and USA breast and plastic reconstructive surgeons. *ecancer*. 2019. DOI: 10.3332/ecancer.2019.927
26. Ballesio L, et al. Postsurgical Ultrasound Evaluation of Patients with Prosthesis in Acellular Dermal Matrix: Results from Monocentric Experience. *International Journal of Surgical Oncology*. 2019. DOI: 10.1155/2019/7437324
27. Khalil HH, et al. Nipple-Sparing Mastectomy and Prepectoral Implant/Acellular Dermal Matrix Wrap Reconstruction in Large Ptotic Breasts. *Plastic and Reconstructive Surgery - Global Open*. 2019. DOI: 10.1097/GOX.0000000000002289

28. Vidya R, et al. Rippling Associated with Pre-Pectoral Implant Based Breast Reconstruction: A New Grading System. *World J Plast Surg.* 2019. DOI: 10.29252/wjps.8.3.311
29. Onesti MG, et al. ADM-assisted prepectoral breast reconstruction and skin reduction mastectomy: Expanding the indications for subcutaneous reconstruction. *Journal of Plastic, Reconstructive & Aesthetic Surgery.* 2019. DOI: 10.1016/j.bjps.2019.11.004
30. Mangialardi ML, et al. Delayed acellular dermal matrix assisted prepectoral breast reconstruction: Preliminary results. *Breast J.* 2019. DOI: 10.1111/tbj.13670
31. Knight HJ, et al. Significantly reducing implant loss rates in immediate implant-based breast reconstruction: A protocol and completed audit of quality assurance. *Journal of Plastic, Reconstructive & Aesthetic Surgery.* 2019. DOI: 10.1016/j.bjps.2019.12.005
32. Chandarana M, Harries S, on behalf of the National Braxon Audit Study Group. Multicentre study of prepectoral breast reconstruction using acellular dermal matrix. *BJS Open.* 2020. DOI: 10.1002/bjs5.50236
33. Cuomo R, et al. Optimization of Prepectoral Breast Reconstruction. *Breast Care.* 2020. DOI: 10.1159/000506347
34. Polotto S, et al. One-step prepectoral breast reconstruction with porcine dermal matrix-covered implant: a protective technique improving the outcome in post-mastectomy radiation therapy setting. *Gland Surg.* 2020. DOI: 10.21037/gs.2020.01.16
35. Maruccia M, et al. Skin-Reducing Mastectomy and Pre-pectoral Breast Reconstruction in Large Ptotic Breasts. *Aesth Plast Surg.* 2020. DOI: 10.1007/s00266-020-01616-2
36. Masiá J & the iBAG Working Group. The largest multicentre data collection on prepectoral breast reconstruction: The iBAG study. *J Surg Oncol.* 2020. DOI: 10.1002/jso.26073
37. Ribuffo D, et al. Dual-Plane Retro-pectoral Versus Pre-pectoral DTI Breast Reconstruction: An Italian Multicenter Experience. *Aesth Plast Surg.* 2020. DOI: 10.1007/s00266-020-01892-y
38. Naemonitou F, et al. Outcome of complete acellular dermal matrix wrap with polyurethane implant in immediate prepectoral breast reconstruction. *Arch Plast Surg.* 2020. DOI: 10.5999/aps.2020.01207
39. Caputo GG, et al. Quality of Life and Early Functional Evaluation in Direct-to-Implant Breast Reconstruction After Mastectomy: A Comparative Study Between Prepectoral Versus Dual-Plane Reconstruction. *Clinical Breast Cancer.* 2020. DOI: 10.1016/j.clbc.2020.11.013
40. Mangialardi ML, et al. Prepectoral implant pocket conversion in breast reconstruction. *JPRAS Open.* 2020. DOI: 10.1016/j.jprra.2020.08.001
41. Caputo GG, et al. Seroma Formation in Pre-pectoral Implant-Based ADM Assisted Breast Reconstruction: A Comprehensive Review of Current Literature. *chr.* 2021. DOI: 10.21614/chirurgia.116.2Suppl.S16
42. Spengler C, et al. Lessons Learned from Three Different Acellular Dermal Matrices in Direct-to-Implant Breast Reconstruction. *MPS.* 2021. DOI: 10.4236/mps.2021.111004
43. Fanelli B, et al. Thoracic migration of silicone gel after breast implant rupture: a case report and literature review. *Eur J Plast Surg.* 2021. DOI: 10.1007/s00238-021-01799-0
44. Mathew J. Can we safely accommodate larger volume implants in inframammary fold nipple sparing mastectomy compared to nipple sacrificing mastectomy in implant-based reconstruction with acellular dermal matrix? *JPRAS Open.* 2021. DOI: 10.1016/j.jprra.2020.10.005
45. Mura S, et al. Direct-to-implant, prepectoral breast reconstruction with Braxon® dermal matrix: A single-center experience with 111 cases. *Breast J.* 2021. DOI: 10.1111/tbj.14193
46. Bojanic C, et al. First use of Braxon® acellular dermal matrix for complex revision aesthetic breast surgery—revision augmentation mastopexy. *Journal of Surgical Case Reports.* 2021. DOI: 10.1093/jscr/rjab256
47. Mathew J. Short- to Medium-term Outcome of Prepectoral versus Subpectoral Direct-to-implant Reconstruction using Acellular Dermal Matrix. *Plastic and Reconstructive Surgery - Global Open.* 2021. DOI: 10.1097/GOX.00000000000003747
48. Maruccia M, et al. Prepectoral breast reconstruction: an ideal approach to bilateral risk-reducing mastectomy. *Gland Surg.* 2021. DOI: 10.21037/gs-21-339
49. Gardani M, et al. Skin-reducing mastectomy and prepectoral breast reconstruction using the Braxon® ADM: a single-centre experience. *Minerva Surg.* 2021. DOI: 10.23736/S2724-5691.21.08941-3
50. Saibene T, et al. Incidence of capsular contracture on irradiated acellular dermal matrices (ADMs)-assisted prepectoral breast reconstructions. *Ann Breast Surg.* 2022. DOI: 10.21037/abs-21-141
51. Quintero Sierra LA, et al. Tissue-Material Integration and Biostimulation Study of Collagen Acellular Matrices. *Tissue Eng Regen Med.* 2022; DOI: 10.1007/s13770-021-00420-6
52. Bassetto F, et al. Complete Implant Wrapping with Porcine-Derived Acellular Dermal Matrix for the Treatment of Capsular Contracture in Breast Reconstruction: A Case–Control Study. *Aesth Plast Surg.* 2022; DOI: 10.1007/s00266-022-02826-6
53. Berna G, et al. Novel three-dimensional acellular dermal matrix for prepectoral breast reconstruction: First year in review with BRAXON®Fast. *Front Surg.* 2022. DOI: 10.3389/fsurg.2022.970053
54. Bassetto F, et al. Braxon®-assisted prepectoral breast reconstruction: A decade later. *Front Surg.* 2022. DOI: 10.3389/fsurg.2022.1009356
55. Gurrado A, et al. Mastectomy with one-stage or two-stage reconstruction in breast cancer: analysis of early outcomes and patient's satisfaction. *Updates Surg.* 2022. DOI: 10.1007/s13304-022-01416-0
56. Caputo GG, et al. From Sub-Pectoral to Pre-Pectoral Implant Reconstruction: A Decisional Algorithm to Optimise Outcomes of Breast Replacement Surgery. *Healthcare.* 2023. DOI: 10.3390/healthcare11050671
57. Santanelli di Pompeo F, et al. Immediate Prepectoral Breast Reconstruction Using an ADM with Smooth Round Implants – A Prospective Observational Cohort Study. *Journal of Plastic, Reconstructive & Aesthetic Surgery.* 2023. DOI: 10.1016/j.bjps.2023.02.014
58. Caputo GG, Mura S, Contessi Negrini F, Albanese R, Parodi PC. From Sub-Pectoral to Pre-Pectoral Implant Reconstruction: A Decisional Algorithm to Optimise Outcomes of Breast Replacement Surgery. *Healthcare.* 2023 Feb 24;11(5):671. doi: 10.3390/healthcare11050671
59. Caputo GG, et al. Optimizing Acellular Dermal Matrix Integration in Heterologous Breast Reconstructive Surgery: Surgical Tips and Post-Operative Management. *Medicina.* 2023 Jun 30;59(7):1231. doi: 10.3390/medicina59071231
60. Cazzato V, Bottosso S, Marchi G, Renzi N, Papa G. A New Ultrasound-assisted Protocol to Avoid Seroma-related Complications in Prepectoral Breast Reconstruction with Acellular Dermal Matrix. *Plastic and Reconstructive Surgery - Global Open.* 2023 Jul;11(7):e5145. doi: 10.1097/GOX.00000000000005145
61. Polotto S, Pedrazzi G, Bergamini M, D'Abbio N, Cattelan L. ADM-assisted Direct-to-implant prepectoral breast reconstruction in post-mastectomy radiation therapy setting: long-term results. *Clinical Breast Cancer.* 2023 Jul;S1526820923001684. doi: 10.1016/j.clbc.2023.06.011

Patented by



**DECOMED**<sup>®</sup>

BIOLOGICAL IDENTITY

info@decomed.it

---

[www.braxon.com](http://www.braxon.com)